

Avaco

Vacuum therapy



Features

product code

manual mode

independent treatment channels

A-UV-AST-VCWH



1

Vacuum therapy

operation modes: continuous, pulse, variable pulse, synchronized with the electric current, pulse with increased lower vacuum	~
two independent treatment channels (easy operation in single-channel mode)	\checkmark
self-sealing suction cups	\checkmark
possibility of vacuum adjustment	\checkmark
current parameter control during the treatment procedure	\checkmark
container filling signalization	\checkmark
dedicated mode for emptying the water container	\checkmark
independent sockets for the application of electrotherapy electrodes	\checkmark

Vacuum therapy technical parameters

vacuum wave pulse frequency	6, 12, 20, 40, 60 pulse/min
vacuum	max 500 mbar/hPa
vacuum adjustment range	100 - 500 mbar, step 40 mbar
number of suction cups	4

General technical parameters

nin	dimensions	30 x 26 x 12 cm
Pa	device weight	3,5 kg
40 bar	power supply, power consumption	230 V, 50 Hz, 20 W, 40 VA

