



Ultrasound head GS-1cm² / 1MHz &3,5 MHz

Ultrasound therapy



Features

product code	A-AS-AST-GS1WH
waterproof ultrasound heads	✓
continuous / pulse emission	✓
ultrasound head contact control (effective treatment time measured)	✓
head sensitivity calibration according to the needs	✓

Ultrasound therapy technical parameters

operating frequency	1 i 3,5 Mhz
effective radiation area	1 cm ²

General technical parameters

accessory dimensions without connection cable	5,5 x 3,7 x 18,4 cm
weight	0,8 kg

Cooperates with

Sonaris S	✓
PHG 200A	✓
PHG 300A	✓
PHG 601C	✓
PHG 701C	✓
PHG 701I	✓