



R 80mW / 660nm red laser point probe

Biostimulation laser therapy



Features

product code	A-AL-AST-80RDV3
operation signalization indicator	✓
optical fiber applicators for laserpuncture and ENT applications	✓

Laser therapy technical parameters

wavelength	660 nm
maximum power	80 mW

General technical parameters

accessory dimensions without connection cable	Ø 1,95 x 16,3 cm
weight	0,5 kg

Cooperates with

Polaris 2	✓
Etius LM	✓
Etius ULM	✓
Polaris HP S	✓
Polaris HP M	✓
PHG 400C	✓
PHG 500I	✓
PHG 601C	✓
PHG 701C	✓
PHG 701I	✓
PhysioGo.Lite LASER	✓