



CL1800 cluster laser applicator

Biostimulation laser therapy



Features

product code	A-AL-AST-CL1800WH
operation signalization indicator	✓
possibility of different operation modes selection : R 5 x 40 mW cluster, 4 x 400 mW IR cluster, 5 x 40 mW R + IR cluster + 4 x 400 mW	✓

Laser therapy technical parameters

wavelength	4x 808 nm and 5x 660 nm
maximum power	4x 400 mW and 5x 40 mW

General technical parameters

accessory dimensions without connection cable	6,2 x 6,8 x 22,5 cm
weight	0.5 kg

Cooperates with

Polaris 2	✓
Etius LM	✓
Etius ULM	✓
Polaris HP S	✓
Polaris HP M	✓
PHG 400C	✓
PHG 500I	✓
PHG 601C	✓
PHG 701C	✓
PHG 701I	✓
PhysioGo.Lite LASER	✓