



Scanning laser applicator PM2-SK2 R-IR 100 / 450mW

Biostimulation laser therapy



Features

product code	A-AL-ASTSK450V2WH
operation signalization indicator	✓
high power of radiation sources measured at the output	✓
three modes of treatment field irradiation in scanning laser applicators	✓
automatic calculation of time relative to treatment parameters - dose, power, duty factor, treatment area	✓
running system with brakes	✓
shelf for the controller and glasses	✓

Laser therapy technical parameters

wavelength	808 nm and 660 nm
maximum power	450 mW and 100 mW

General technical parameters

accessory dimensions without connection cable	15,0 x 7,0 x 12,2 cm
size area needed for the stand base	69,0 x 69,0 cm
arm height adjustment	60,0 - 140,0 cm
scanning laser applicator rotation adjustment angle	-90° - +90°
weight	1 kg
weight with a stand	12,5 kg

Cooperates with

Polaris 2	✓
Etius LM	✓
Etius ULM	✓
Polaris HP S	✓
Polaris HP M	✓
PHG 400C	✓
PHG 500I	✓
PHG 601C	✓
PHG 701C	✓
PHG 701I	✓
PhysioGo.Lite LASER	✓