



Plate magnetic field applicator - CPE2 type

Magnetotherapy



Features

product code	A-AG-AST-CPE2WH
operation with CPE1 plate applicator (CPE2 does not work alone)	✓
convenient application of applicators with straps and velcro belts	✓

Magnetotherapy technical parameters

maximum magnetic field induction	10 mT
----------------------------------	-------

General technical parameters

accessory dimensions without connection cable	10 x 5,5 x 15 cm
weight	1,2 kg

Cooperates with

Etius LM	✓
Etius ULM	✓
PHG 500I	✓
PHG 701I	✓