



CS75 solenoid applicator

Magnetotherapy



Features

product code	A-AG-AST-CS75A
magnetic field activity visual indicator in the form of an illuminator	✓
magnetic field activity sensory indicators in the form of magnetic stripes	✓

Magnetotherapy technical parameters

maximum induction in the geometrical center of the applicator:	3 mT
maximum induction at the applicator site	4,5 mT
maximum induction change value	9 mT

General technical parameters

accessory dimensions without connection cable	Ø 83,2 x 24,2 cm
weight	18 kg (Al) / 25 kg (Cu)

Cooperates with

PMG 827	✓
---------	---