



CS60 solenoid applicator

Magnetotherapy



Features

product code	A-AG-AST-CS60A
magnetic field activity visual indicator in the form of an illuminator	✓
magnetic field activity sensory indicators in the form of magnetic stripes	✓

Magnetotherapy technical parameters

maximum induction in the geometrical center of the applicator:	4 mT
maximum induction at the applicator site	6 mT
maximum induction change value	12 mT

General technical parameters

accessory dimensions without connection cable	Ø 68,2 x 24,2 cm
weight	16 kg (Al) / 22 kg (Cu)

Cooperates with

PMG 815	✓
PMG 825	✓
PMG 827	✓