



CPEP plate applicator

Magnetotherapy



Features

product code	A-AG-AST-CPEP	
convenient application of applicators with straps and velcro belts		✓
magnetic field activity sensory indicators in the form of magnetic stripes		✓

Magnetotherapy technical parameters

maximum induction in the geometrical center of the applicator:	25 mT
maximum induction at the applicator site	50 mT
maximum induction change value	100 mT

General technical parameters

accessory dimensions without connection cable	10,0 x 5,5 x 15,0 cm
weight	2,2 kg

Cooperates with

PMG 825	✓
PMG 827	✓