



# CP plate applicator

Magnetotherapy



## Features

product code	A-AG-AST-CPV2
magnetic field activity visual indicator in the form of an illuminator	✓
magnetic field activity sensory indicators in the form of magnetic stripes	✓
running system with brakes	✓

## Magnetotherapy technical parameters

maximum induction in the geometrical center of the applicator:	3 mT
maximum induction at the applicator site	12,5 mT
maximum induction change value	25 mT

## General technical parameters

accessory dimensions without connection cable	20,0 x 6,9 x 15,0 cm
size area needed for the stand base	65,0 x 75,0 cm
stand height adjustment	40,0 - 110,0 cm
applicator rotation adjustment angle	150°
weight	5,5 kg
weight with a stand	29 kg

## Cooperates with

PMG 815	✓
PMG 825	✓
PMG 827	✓