



# Versa trolley

## Trolleys

### Features

|                                                                    |              |
|--------------------------------------------------------------------|--------------|
| product code                                                       | A-AM-AST-VSA |
| compatibility with ASTAR units / accessories                       | ✓            |
| running system with brakes                                         | ✓            |
| shelves height adjustment (second, third and for CPEP applicators) | ✓            |

### General technical parameters

|                                                                    |                    |
|--------------------------------------------------------------------|--------------------|
| upper shelf load capacity                                          | max 10 kg          |
| total load capacity of the bottom shelf (with drawer and contents) | max 7 kg           |
| bottom shelf adjustment range                                      | approx. 60 cm      |
| external dimension of the upper shelf (WxD)                        | 58 x 34 cm         |
| internal dimension of the upper shelf (WxD)                        | 39,5 x 30 cm       |
| external dimension of the bottom shelf (WxD)                       | 40 x 26 cm         |
| external dimension of the drawer (WxDxH)                           | 39 x 33 x 9 cm     |
| internal size of the drawer (WxDxH)                                | 28,5 x 18,5 x 8 cm |
| trolley size (WxDxH)                                               | 58 x 49 x 87,5 cm  |
| weight                                                             | 13,9 kg            |

### Dedicated for

|              |   |
|--------------|---|
| Sonaris M    | ✓ |
| Sonaris S    | ✓ |
| Polaris 2    | ✓ |
| Etius        | ✓ |
| Etius U      | ✓ |
| Etius LM     | ✓ |
| Etius ULM    | ✓ |
| Polaris HP S | ✓ |
| Polaris HP M | ✓ |
| PHG 100A     | ✓ |

|                       |   |
|-----------------------|---|
| PHG 200A              | ✓ |
| PHG 300A              | ✓ |
| PHG 400C              | ✓ |
| PHG 500I              | ✓ |
| PHG 601C              | ✓ |
| PHG 701C              | ✓ |
| PHG 701I              | ✓ |
| PhysioGo.Lite LASER   | ✓ |
| PhysioGo.Lite SONO    | ✓ |
| PhysioGo.Lite ELECTRO | ✓ |