



Ultrasound head GU-5cm2 / 1MHz & 3 MHz

Ultrasound therapy



Features

product code	A-AS-AST-GU5
LIPUS mode	✓
water resistance (IPX7)	✓
continuous / pulse emission	✓
lack of contact detection	✓
US head temperature control	✓
US head sensitivity adjustment	✓

Ultrasound therapy technical parameters

acoustic working frequency	1 MHz, 3 MHz
acoustic working frequency - switch mode	1/3 MHz (switching frequency every 8 sec.)
acoustic working frequency - LIPUS mode	1 MHz
total area of the head front	5 cm ²
effective radiation area	3,4 cm ²

General technical parameters

accessory dimensions without connection cable	5,5 x 3,7 x 18,4 cm
weight	max 0,5 kg

Cooperates with

PhysioGo.Lite SONO	✓
PhysioGo.Lite COMBO	✓