



Carro Versa X

Carros

Característica

código de producto

compatibilidad con los dispositivos / accesorios ASTAR

sistema de tracción con frenos

regulación de la altura de las repisas (segunda, tercera y para los aplicadores CPEP)

A-AM-AST-VSX



Parámetros técnicos generales

capacidad de carga de la repisa de superior	máx 10 kg
capacidad de carga total de la repisa inferior (con el cajón y el contenido)	máx 7 kg
rango de regulación de la repisa inferior	approx 50 cm
dimensiones externas de la repisa superior (AxL)	58,0 x 39,0 cm
dimensiones internas de la repisa superior (AxL)	39,5 x 35,0 cm
dimensiones externas de la repisa inferior (AxL)	40,0 x 26,0 cm
dimensiones externas del cajón (AxLxL)	39,0 x 33,0 x 9,0 cm
dimensiones internas del cajón (AxLxL)	28,5 x 18,5 x 8,0 cm
dimensiones del carro (AxLxA)	58,0 x 49,0 x 89,0 cm
peso	23 kg

Dedicado a

Sonaris M	✓
Sonaris S	✓
Polaris 2	✓
Etius	✓
Etius U	✓
Etius LM	✓
Etius ULM	✓
Polaris HP S	✓
Polaris HP M	✓
Impactis M	✓

Impactis M+	✓
PHG 100A	✓
PHG 200A	✓
PHG 300A	✓
PHG 400C	✓
PHG 500I	✓
PHG 601C	✓
PHG 701C	✓
PHG 701I	✓
PhysioGo.Lite LASER	✓
PhysioGo.Lite SONO	✓
PhysioGo.Lite ELECTRO	✓